

# simpleABI<sup>TM</sup>

Available with optional advanced vascular testing!

## simpleABI<sup>TM</sup> Systems

*Accurate, Low Cost System for The Diagnosis and Documentation of P.A.D.*

The fast, easy way to create the reports needed for ABI and advanced vascular examinations. No secondary transferring of data or hidden costs of ownership.\*



**Newman Medical** LOWER EXTREMITY PHYSIOLOGIC STUDY SINGLE LEVEL  
ABI EXAM

Patient Name: Paul Smith, Patient ID: AN22956  
 Patient Age: 69, Date: 11/20/11  
 Physician: Dr. Brody, Examiner: Alice Johnson

**Risk Factors**  
 Tobacco Use,  Hypertension,  Diabetes,  Heart Disease,  Prev CV Event,  Hyperlipidemia,  Stroke/TIA,  Prev Vasc Surgery,  Cholesterol > 200,  Cholesterol > 240

**Current Symptoms**  
 Intermittent Claudication,  Numbness/Tingling in Feet,  Leg Cramps,  Leg Pain,  Gangrene,  Location: Subtoses,  Thigh,  Calf,  Heel

**ABI Guidelines**  
 >1.30: Incompressible  
 1.00-1.29: Normal  
 0.91-0.99: Borderline  
 0.41-0.90: Moderate/Severe  
 0.00-0.40: Severe

**Comments:**

**RIGHT**  
 BRACHIAL: 124, 138  
 RIGHT ABI: 1.13

**LEFT**  
 BRACHIAL: 118, 142  
 LEFT ABI: 1.16

**CPT CODE:**  
 93921 Athroangioplasty of the artery with intermittent claudication  
 93920 Diabetes Mellitus with Peripheral Circulatory Disorders  
 93922 Non-invasive physiologic studies of the upper or lower extremity arteries, single level, bilateral

www.newman-medical.com 800-267-5549

The simpleABI Systems from Newman Medical use the power of the newest portable computers to produce the reimbursement reports you need quickly, easily, and accurately. The diagnosis is made using the gold standard Doppler and simple pulse volume recording (PVR) waveforms. Working directly on the computer is simple, intuitive, effective, and reduces the time needed to get to the final report necessary for reimbursement purposes. The system meets the new requirements\* for CPT code 93922, 93923 & 93924.

\*New CPT codes for 2011 now require the measurement of the ABI to be done with BOTH the posterior tibial and the dorsalis pedis arteries. Many systems are not capable of this and may no longer be eligible for reimbursement.

Exams/Week	Avg. Reimbursement	Monthly Revenue	Annual Revenue
1	\$114	\$494	\$5,923
5	\$114	\$2,468	\$29,617
10	\$114	\$4,936	\$59,234



**Newman Medical<sup>TM</sup>**

Call us at **800.267.5549**

[www.newman-medical.com](http://www.newman-medical.com)

42 Sherwood Terrace, Suite 2 | Lake Bluff, IL 60044

## simpleABI<sup>®</sup> Reimbursable Vascular Systems

### Complete Systems for fast, accurate, reimbursable vascular diagnosis

- A compact and powerful computer and wireless printer for reimbursable report generation
- A DigiDop model DD-301 Doppler with a sensitive, comfortable DD-8MHz vascular probe  
*(5MHz also available)*
- 4 special vascular cuffs and aneroid  
*(4 additional leg cuffs, toe cuff, & PPG with simpleABI+)*
- Sturdy rolling stand with basket
- Long life rechargeable batteries  
Up to 10 hours of continuous operation on battery power
- Pulse Volume Recording Hardware for reimbursement purposes
- 2 patient informational posters for office exam rooms
- 200 Patient flyers with risk assessment quiz
- Operating manual and ABI instructions
- American Heart Association P.A.D. Guidelines



- **ABI-300** | simpleABI System - PC based ABI system - Includes: PC, wireless printer, 8MHz Doppler, Dual PVR, 4 cuffs, aneroid and roll stand
- **ABI-300-NOPC** | simpleABI without computer, printer and stand
- **ABI-500** | simpleABI+ System - PC based system for advanced vascular testing. Includes: ABI-300, simpleABI+ software upgrade, 4 additional leg cuffs, toe cuff, PPG
- **ABI-500-NOPC** | simpleABI+, without computer, printer and stand
- **ABI-501** | simpleABI+ field upgrade for previously purchased ABI-300 Systems

Call us at **800.267.5549** or Fax to **847.283.9189**

visit us online at [www.newman-medical.com](http://www.newman-medical.com)

Made in the U.S.A!

NM0022-C