

# simpleABI<sup>TM</sup>

## simpleABI<sup>TM</sup> System 300

*Accurate, Low Cost System for  
The Diagnosis and Documentation  
of P.A.D.*

The fast, easy way to create the reports needed for ABI examinations. No secondary transferring of data or hidden costs of ownership.\*

**NEWMAN Medical** LOWER EXTREMITY PHYSIOLOGIC STUDY SINGLE LEVEL  
ABI EXAM

Patient Name: Paul Smith, Patient ID: A022946  
 Patient Age: 62, Date: 11/20/11  
 Physician: Dr. Brodie, Examiner: Alice Johnson

**Risk Factors:**  
 Tobacco Use,  Hypertension,  Diabetes,  Heart Disease,  Prev CV Event,  Hyperlipidemia,  Stroke/TIA,  Prev Vasc Surgery,  Cholesterol > 70

**Current Symptoms:**  
 Intermittent Claudication,  Numbness/Tingling in Feet,  Leg Pain,  Gait Issues,  Location: Subtoses,  Thigh,  Calf,  Heel

**ABI Guidelines:**  
 >1.30: Incompressible  
 1.00-1.29: Normal  
 0.91-0.99: Borderline  
 0.41-0.90: Moderate/Severe  
 0.00-0.40: Severe

**ABI Results:**  
 RIGHT BRACHIAL: 124, LEFT BRACHIAL: 118  
 RIGHT DP: 140, LEFT DP: 144  
 RIGHT ABI: 1.13, LEFT ABI: 1.16

**ICD-9 CODE:**  
 440.21 Atherosclerosis of the artery with intermittent claudication  
 250.70 Diabetes Mellitus with Peripheral Circulatory Disorders  
 441.90 Peripheral vascular disease, unspecified

Comments: \_\_\_\_\_

CPT CODE: 93922 Non-invasive physiologic studies of the upper or lower extremity arteries, single level, bilateral.

www.newman-medical.com 800-267-5549



The simpleABI System 300 from Newman Medical uses the power of the newest portable computers to produce the reimbursement reports you need quickly, easily, and accurately. The diagnosis is made using the gold standard Doppler ABI and simple pulse volume recording (PVR) waveforms. Working directly on the computer is simple, intuitive, effective, and reduces the time needed to get to the final report necessary for reimbursement purposes. The system meets the new requirements\* for CPT code 93922.

\*New CPT codes for 2011 now require the measurement of the ABI to be done with BOTH the posterior tibial and the dorsalis pedis arteries. Many systems are not capable of this and may no longer be eligible for reimbursement.

Exams/Week	Avg. Reimbursement	Monthly Revenue	Annual Revenue
1	\$114	\$494	\$5,923
5	\$114	\$2,468	\$29,617
10	\$114	\$4,936	\$59,234

\*Competitive systems may cost you an additional \$800 PER YEAR with as few as 5 exams per week – see [www.newman-medical.com](http://www.newman-medical.com) for details.



**Newman Medical<sup>TM</sup>**

Call us at **800.267.5549**

[www.newman-medical.com](http://www.newman-medical.com)

42 Sherwood Terrace, Suite 2 | Lake Bluff, IL 60044

## simpleABI Reimbursable Vascular System 300

### COMPLETE SYSTEM FOR FAST, ACCURATE, REIMBURSABLE DIAGNOSIS OF P.A.D

- A compact and powerful HP netbook and wireless printer for reimbursable report generation.
- A DigiDop model DD-301 Doppler with a sensitive, comfortable DD-8MHz vascular probe
- 4 special vascular cuffs and aneroid.
- Sturdy rolling stand with basket
- Long life rechargeable batteries. Up to 10 hours of continuous operation on battery power.
- Pulse Volume Recording Hardware for reimbursement purposes.
- 2 patient informational posters for office exam rooms
- 200 Patient flyers with risk assessment quiz
- Operating manual and ABI instructions
- American Heart Association P.A.D. Guidelines
- If you have a laptop, ask about our System 100.



Call us at **800.267.5549** or Fax to **847.283.9189**

visit us online at [www.newman-medical.com](http://www.newman-medical.com)

Made in the U.S.A!

NM0022-B